PULSE OXIMETER

EN - INSTRUCTION FOR USE

INTRODUCTION
Thank you for choosing the LANAFORM™ "PULSE OXIMETER S1". This Pulse Oximeter is an easy, reliable, non-invasive and continuous way of measuring your blood oxygen saturation and pulse rate, whether you're an adult or a child.

A FEW WORDS ABOUT THE PULSE OXIMETER...

1- **Uses of this device**
The Pulse Oximeter can be used to measure your functional blood haemoglobin oxygen saturation (SpO2) and heart rate. It works when placed on your finger and can be used at home as well in a sporting context by socio-medical organisations, etc. This device is not intended to be used under continuous supervision.

2- **General Description**
The Pulse Oximeter is a non-invasive method used to blood oxygen saturation (SpO2) by checking the percentage of oxygen-saturated haemoglobin and measuring heart rate. This process has been in regular use now for twenty years and is proving to be an essential control in medical practice for maintaining oxygen levels and preventing respiratory problems. The value of the oxygen saturation measured is regarded as dangerous if it is below 90% for 90% to 100% under normal anatomical conditions.

The advantage of the LANAFORM oximeter is that takes up little space, has low energy consumption, is easy to handle and carry. By placing your finger in the photoelectric sensor, the analysis of the haemoglobin-saturated values measured and the heart rate are displayed immediately.

3- **Measurement principle**
The measurement principle is based on the de Beer-Lambert Act. The principle is based on two lights emitted (red and infrared) of 660 and 940nm respectively and of the measurement of their absorption by the pulsatile flow. The sensor consists of two diodes emitting a red light, which must be facing an easily identifiable receptive field. The best results are obtained by placing the transmitter on the nail. The information acquired will thus be displayed by two series of LEDs, which process information through two microprocessors.

TITLE AND DESCRIPTION OF DISPLAY SYMBOLS ON THE BROCHURE DESCRIBING THE PRODUCT'S FEATURES:

Symbol	Definition	Log in Register	Information provided by the manufacturer
Manufacturer's details			LANAFORM SA B-4141 SPRIMONT
Batch number		LOT	
User Guide			Follow the instructions for use
BF device			Please note: consult the enclosed documents. Please read these instructions before using the device. Please also follow these instructions carefully when using the device.
Please note: see instructions for use			Please note: consult the enclosed documents. Please read these instructions before using the device. Please also follow these instructions carefully when using the device.
Notified body			
Electrical and electronic waste (WEEE)			
Device not intended for continuous monitoring (no SpO2 alarm)			
Oxygen saturation percentage display		%SpO ₂	
Pulse rate display			

PLEASE READ ALL THE INSTRUCTIONS BEFORE USING YOUR "PULSE OXIMETER S1", ESPECIALLY THIS BASIC SAFETY INFORMATION:

- Do not use the pulse oximeter in an MRI room.
- Keep the oximeter away from young children. The battery compartment cover and the strap can cause strangulation or choking.
- Some activities will expose you to injury, including the risk of strangulation, if the strap gets wrapped around the neck. Use the strap with care.
- The measurements acquired by this pulse oximeter do not replace medical supervision. It is only a source of extra information to be passed on to your care provider.
- The oximeter may mistakenly interpret excessive movement as a good-quality pulsatile signal. Please stay still when using this device.
- The oximeter must be able to measure the pulse correctly to give an exact reading. Put the device on the same hand/arm as a blood pressure cuff or blood pressure monitoring device.
- Do not put the oximeter in any liquid and do not clean it with cleaners containing ammonium chloride, alcohol or products not mentioned in this user's guide.
- The LANAFORM pulse oximeter is not intended to be used in a medical establishment.
- The oximeter's performance might be impaired in any of the following cases:
 - fluctuating or very bright light
 - weak pulse
 - low haemoglobin rate
 - arterial catheter
 - nail polish and/or false nails
 - recent tests requiring the injection of an intravenous drug.
- The oximeter might not run if you have poor blood circulation. Rub your finger to increase circulation or put the device on another finger.
- The batteries may over-charge or explode if not used correctly or are disposed of in a manner that does not comply with the regulations. Remove the battery if you are intending to store the device for more than 30 days.
- Do not use this device in a room in which aerosols (sprays) are used or in a room in which oxygen is administered.
- Do not use the oximeter outside the specified operating temperature ranges.
- Do not use the device for more than 30 minutes without changing fingers.
- The oximeter must be used according to the instructions in the user manual.
- Do not use the device in the vicinity of large electromagnetic fields, such as cordless or portable telephones.
- Comply with the local recycling laws when disposing of the oximeter and its components, including the batteries.
- Always keep the oximeter away from heat sources such as stoves, radiators, etc.
- If there is a problem with your device, please contact your reseller.
- Do not try to repair the device yourself.
- This equipment is not intended for use by people, including children, with limited physical, sensory or mental abilities, or by people lacking experience or knowledge, unless they are supervised or have been instructed beforehand on how to use the equipment by someone responsible for their safety. Children should be supervised to ensure that they do not play with the device.
- This device must be used at an ambient temperature between 5°C and 40°C.
- Do not expose the device to electric shocks.
- Do not expose this oximeter to extremes of temperature above 60°C or below -20°C.
- Do not use the device if the relative humidity is above 80% or below 30%.

TROUBLE SHOOTING

If you encounter any of the following problems when using your PULSE OXIMETER S1, please consult this troubleshooting guide to resolve the problem. If the problem persists, please contact your reseller.

Anomalies	Possible causes	Solutions
SpO2 and pulse rate display is unstable	The cuff is not positioned correctly inside the device. Your finger moves or the patient is moving. The batteries are not charged.	Put your finger in the correct place and try again. Stay still during the measuring process. Change the batteries.
The device does not light up	The batteries are not inserted correctly.	Re-insert the batteries.
The device is malfunctioning	The device is malfunctioning. catheter arterial - catheter arterial - venous - doigts à ongles et/ou faux ongles	Contact the supplier or its after-sales department. Normal
The light suddenly goes out.	The batteries are nearly flat.	Change the batteries.

WARRANTY

EC 0029 DECLARATION OF CONFORMITY

We, the company LANAFORM, hereby declare that the product LANAFORM Batch 001 medical device complies with European Directive 93/42/EC concerning medical devices and its most recent amendment 2007/47/EC. APRAGAZ 156 chaussees de Wivoroode 1120 Brussels (Belgium), notified body No. 0029, has issued the Certificate of Conformity pursuant to Appendix IV in...
We, the company LANAFORM, hereby declare that the product LANAFORM Batch 001 medical device complies with European Directive 93/42/EC concerning medical devices and its most recent amendment 2007/47/EC. APRAGAZ 156 chaussees de Wivoroode 1120 Brussels (Belgium), notified body No. 0029, has issued the Certificate of Conformity pursuant to Appendix IV in...

ADVICE REGARDING THE DISPOSAL OF WASTE

- All the packaging is composed of materials that pose no hazard for the environment and which can be disposed of at your local sorting centre to be used as secondary raw materials. The cardboard may be disposed of in a paper recycling bin. The packaging film must be taken to your local sorting and recycling centre.
- When you have finished using the device, please dispose of it in an environmentally friendly way and in accordance with the law. Prior to disposal, please remove the batteries and dispose of it at a collection point so that it can be recycled.

Please note: under no circumstances should spent batteries be disposed of with household waste!

LIMITED WARRANTY

LANAFORM™ guarantees this product against any material or manufacturing defect for a period of two years from the date of purchase, except in the circumstances described below.
The LANAFORM™ guarantee does not cover damage caused as a result of normal wear to this product. In addition, the guarantee covering this LANAFORM™ product does not cover damage caused by abusive or inappropriate or incorrect use, accidents, the use of unauthorized accessories, changes made to the product or any other circumstance, of whatever sort, that is outside LANAFORM™'s control.
LANAFORM™ may not be held liable for any type of circumstantial, indirect or specific damage.
All implied warranties relating to the suitability of the product are limited to a period of two years from the initial date of purchase as long as a copy of the proof purchase can be supplied.
On receipt, LANAFORM™ will repair or replace your appliance at its discretion and will return it to you. The guarantee is only valid through LANAFORM™'s Service Centre*. Any attempt to maintain this product by a person other than LANAFORM™'s Service Centre will render this guarantee void.

DESCRIPTION OF THE OXIMETER DISPLAY

The display interface of the oximeter can automatically point in four directions the display of data automatically according to the parameters detected. There are four types of possible display all based on the four directions.

- FUNCTIONS**
- 1. LED - colour display, four display modes.
 - 2. Automatic display function.
 - 3. Real-time visual alarm function.
 - 4. Low power consumption. Lasts for 50 hours continuously.
 - 5. Low power consumption.
 - 6. Low power consumption.
 - 7. Low power consumption.
 - 8. Low power consumption.
 - 9. Low power consumption.
 - 10. Low power consumption.

TECHNICAL SPECIFICATIONS

Normal conditions of use	
Display resolution	
-SpO2	+/- 1%
-Heart rate	+/- 1 BPM
Operating temperature	0 to 40°C (Moisture content: 15-1<95%)
Storage temperature	-20C-60C (Humidity content: 5%-95%)
Power supply	3V continuous (2 x AAA batteries, not included)
Dimensions	64.5x37.5x35mm (LxWxH)
Weight	Net 36g

Measurement range	
SpO2	70%-99%
Heart rate	30-240BPM


Precision	
SpO2	±2% (70%-99%), not specified (≤70%)
Heart rate	±1 bpm

Default alarm value	
Upper limit: 240BPM	
Lower limit: 30BPM	

Low perfusion	
≤0.4%	
Switches off automatically	8 sec.

INSTRUCTIONS FOR USE

Insert two AAA batteries (not included) into the case, checking that the polarities are correct before replacing the cover.

MAINTENANCE AND STORAGE

- Use an alcohol solution for medical use to clean the rubber sensor in contact with your finger inside the oximeter and clean this section before and after each use.
- Gently dry the oximeter with a clean, dry cloth. Do not use alcohol or other solvents.
- For optimum results, please use the device in a clean, dry environment.
- Please clean your OLED screen using a slightly damp cloth.
- Please remove the battery when not using your device for an long length of time.

PLEASE READ ALL THE INSTRUCTIONS BEFORE USING YOUR "PULSE OXIMETER S1", ESPECIALLY THIS BASIC SAFETY INFORMATION:

1. Do not use the pulse oximeter in an MRI room.

2. Keep the oximeter away from young children. The battery compartment cover and the strap can cause strangulation or choking.

3. Some activities will expose you to injury, including the risk of strangulation, if the strap gets wrapped around the neck. Use the strap with care.

4. The measurements acquired by this pulse oximeter do not replace medical supervision. It is only a source of extra information to be passed on to your care provider.

5. The oximeter may mistakenly interpret excessive movement as a good-quality pulsatile signal. Please stay still when using this device.

6. The oximeter must be able to measure the pulse correctly to give an exact reading. Put the device on the same hand/arm as a blood pressure cuff or blood pressure monitoring device.

7. Do not put the oximeter in any liquid and do not clean it with cleaners containing ammonium chloride, alcohol or products not mentioned in this user's guide.

8. The LANAFORM pulse oximeter is not intended to be used in a medical establishment.

9. The oximeter's performance might be impaired in any of the following cases:

- fluctuating or very bright light
- weak pulse
- low haemoglobin rate
- arterial catheter
- nail polish and/or false nails
- recent tests requiring the injection of an intravenous drug.

10. The oximeter might not run if you have poor blood circulation. Rub your finger to increase circulation or put the device on another finger.

11. The batteries may over-charge or explode if not used correctly or are disposed of in a manner that does not comply with the regulations. Remove the battery if you are intending to store the device for more than 30 days.

12. Do not use this device in a room in which aerosols (sprays) are used or in a room in which oxygen is administered.

13. Do not use the oximeter outside the specified operating temperature ranges.

14. Do not use the device for more than 30 minutes without changing fingers.

15. The oximeter must be used according to the instructions in the user manual.

16. Do not use the device in the vicinity of large electromagnetic fields, such as cordless or portable telephones.

17. Comply with the local recycling laws when disposing of the oximeter and its components, including the batteries.

18. Always keep the oximeter away from heat sources such as stoves, radiators, etc.

19. If there is a problem with your device, please contact your reseller.

20. Do not try to repair the device yourself.

21. This equipment is not intended for use by people, including children, with limited physical, sensory or mental abilities, or by people lacking experience or knowledge, unless they are supervised or have been instructed beforehand on how to use the equipment by someone responsible for their safety. Children should be supervised to ensure that they do not play with the device.

22. This device must be used at an ambient temperature between 5°C and 40°C.

23. Do not expose the device to electric shocks.

24. Do not expose this oximeter to extremes of temperature above 60°C or below -20°C.

25. Do not use the device if the relative humidity is above 80% or below 30%.

Anomalies	Possible causes	Solutions
SpO2 and pulse rate display is unstable	The cuff is not positioned correctly inside the device. Your finger moves or the patient is moving. The batteries are not charged.	Put your finger in the correct place and try again. Stay still during the measuring process. Change the batteries.
The device does not light up	The batteries are not inserted correctly.	Re-insert the batteries.
The device is malfunctioning	The device is malfunctioning. catheter arterial - catheter arterial - venous - doigts à ongles et/ou faux ongles	Contact the supplier or its after-sales department. Normal
The light suddenly goes out.	The batteries are nearly flat.	Change the batteries.

WARRANTY

EC 0029 DECLARATION OF CONFORMITY

We, the company LANAFORM, hereby declare that the product LANAFORM Batch 001 medical device complies with European Directive 93/42/EC concerning medical devices and its most recent amendment 2007/47/EC. APRAGAZ 156 chaussees de Wivoroode 1120 Brussels (Belgium), notified body No. 0029, has issued the Certificate of Conformity pursuant to Appendix IV in...
We, the company LANAFORM, hereby declare that the product LANAFORM Batch 001 medical device complies with European Directive 93/42/EC concerning medical devices and its most recent amendment 2007/47/EC. APRAGAZ 156 chaussees de Wivoroode 1120 Brussels (Belgium), notified body No. 0029, has issued the Certificate of Conformity pursuant to Appendix IV in...

ADVICE REGARDING THE DISPOSAL OF WASTE

- All the packaging is composed of materials that pose no hazard for the environment and which can be disposed of at your local sorting centre to be used as secondary raw materials. The cardboard may be disposed of in a paper recycling bin. The packaging film must be taken to your local sorting and recycling centre.
- When you have finished using the device, please dispose of it in an environmentally friendly way and in accordance with the law. Prior to disposal, please remove the batteries and dispose of it at a collection point so that it can be recycled.

Please note: under no circumstances should spent batteries be disposed of with household waste!

LIMITED WARRANTY

LANAFORM™ guarantees this product against any material or manufacturing defect for a period of two years from the date of purchase, except in the circumstances described below.
The LANAFORM™ guarantee does not cover damage caused as a result of normal wear to this product. In addition, the guarantee covering this LANAFORM™ product does not cover damage caused by abusive or inappropriate or incorrect use, accidents, the use of unauthorized accessories, changes made to the product or any other circumstance, of whatever sort, that is outside LANAFORM™'s control.
LANAFORM™ may not be held liable for any type of circumstantial, indirect or specific damage.
All implied warranties relating to the suitability of the product are limited to a period of two years from the initial date of purchase as long as a copy of the proof purchase can be supplied.
On receipt, LANAFORM™ will repair or replace your appliance at its discretion and will return it to you. The guarantee is only valid through LANAFORM™'s Service Centre*. Any attempt to maintain this product by a person other than LANAFORM™'s Service Centre will render this guarantee void.

DESCRIPTION OF THE OXIMETER DISPLAY

The display interface of the oximeter can automatically point in four directions the display of data automatically according to the parameters detected. There are four types of possible display all based on the four directions.

- FUNCTIONS**
- 1. LED - colour display, four display modes.
 - 2. Automatic display function.
 - 3. Real-time visual alarm function.
 - 4. Low power consumption. Lasts for 50 hours continuously.
 - 5. Low power consumption.
 - 6. Low power consumption.
 - 7. Low power consumption.
 - 8. Low power consumption.
 - 9. Low power consumption.
 - 10. Low power consumption.

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